

Simply Touching!

Enhanced Usage With Surprising Ease-Of-Use!

Tranquility Touch Screen Patient Monitors



T - LITE


8" color TFT display

- ▶ 6 Waveforms displayed
- ▶ Simultaneous multi-lead ECG monitoring
- ▶ Graphical and tabular trending
- ▶ Audible and visual alarms
- ▶ Advanced ST and arrhythmia detection
- ▶ T-Lite Options: Thermal Recorder

TRANQUILITY II

12.1" color TFT display

- ▶ Multi-Lead ECG monitoring
- ▶ Drug dose calculations
- ▶ 4, 6 & 8 waveforms shown
- ▶ NIBP list & OxyCRG dynamic view display
- ▶ Networkable with Central Monitoring System via a wireless or wired network
- ▶ Compact design is convenient for mobile monitoring
- ▶ Optional: Dual IBP, EtCo₂, Anesthetic Agent and Built-In Thermal Recorder



SCHILLER

A M E R I C A